# alcoSense<sup>®</sup> excel

## Manual

Page 1

## Instrukcja

Strona 29

### PROBLEM?

We can solve 99% of issues  
in a quick call 0800 195 0088  
or email [info@alcosense.com](mailto:info@alcosense.com)

### PYTANIA?

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**WARNING:** AlcoSense is a sensitive measurement device. Performing a breath test less than 10 minutes after drinking risks blowing high levels of alcohol directly from your mouth into the system which could permanently damage the sensor and cause the system to read inaccurately. Blood alcohol levels can continue to rise for up to 90 minutes after drinking, as a result we only recommend using the Breathalyser after this time (and ideally the morning after) to ensure your reading will not increase after you test yourself. Do not blow directly into the blow hole without a blow tube – this can damage the sensor and significantly reduce the lifetime of your unit.

# AlcoSense Excel Instruction Manual

## Thank you for purchasing your new AlcoSense Breathalyser.

As you know, you never can tell how long it takes for alcohol to clear from your body after having a drink. It depends on many variables including body mass, time of day, food intake and even if you have a cold.

Road safety is our primary concern. AlcoSense breathalysers are designed to quickly and accurately indicate your level of intoxication, making it easy for you to tell when the alcohol has cleared from your system the morning after a few drinks.

Thank you again for purchasing one of our products and helping to make our roads a safer place. I'm sure you will be delighted with it.

Yours sincerely,

Hunter Abbott, Managing Director

### Important - Please read

It is safest to drive when you have no alcohol left in your blood stream. As a result we only recommend driving when the AlcoSense breathalyser shows zero and you feel 100%. If you have any doubt that you may still have alcohol in your system - do not drive. Before using the AlcoSense breathalyser you agree you have read the operating instructions & are aware of the legal alcohol limit in your jurisdiction. The manufacturer distributors and sellers of this product exclude all liability for improper use or loss or damage caused by actions taken by you after use of this product.

AlcoSense Excel User Manual Version 2.0

## Contents

How Breathalysers Work	4
How Alcohol is Processed by the Body	4
Do's and Don't's	5
Getting Started	6
Unpacking the Box	6
Understanding the Buttons and Screens	7
What the Icons Mean	8
Installing the Batteries	10
First Time Setup	11
How to Perform a Breath Test	12
Understanding the Results	14-15
Changing Settings & Recalling Results	16
Changing the Language	17
Changing Limits, Units of Measure and Driving Abroad	17-18
How can a Breathalyser Display a Blood Alcohol Reading?	19
The Blood Breath Ratio	19
Common Drink Drive Limits and Blood Breath Ratios	20-21
Maintaining your AlcoSense Excel Breathalyser	22
Technical Support	23
Troubleshooting	23-25
Specification	26
Accessories	26
Legal & Warranty	27

## How Breathalysers Work

Breathalysers work by measuring the concentration of alcohol vapour in your breath. Contrary to popular belief the alcohol is not coming from residue in your mouth but is actually evaporating from the blood that passes through your lungs when you breathe. The amount of alcohol evaporating from your lungs is directly proportional to the amount of alcohol in your blood stream making it possible to accurately calculate the per mil Blood Alcohol Concentration (%BAC) or Milligrams of alcohol per Litre of Breath (mg/l) and therefore your level of intoxication.

## How Alcohol Is Processed by the Body

When you drink, alcohol is absorbed into the blood stream through your mouth and throat but mostly through the stomach and small intestine. The more you drink the more alcohol is absorbed by the body and the higher your blood and breath alcohol readings will become. Approximately 90 minutes after your last drink the alcohol is distributed equally around the body and through the blood stream. This is known as the point of Equilibrium.

As blood flows round the body around 90% of the alcohol is broken down in to water and carbon dioxide by the liver in a process called Oxidation. The remaining 10% of the alcohol continues to the lungs and kidneys. A healthy liver typically breaks down alcohol at a rate of one unit of alcohol per hour but sometimes it can take a longer or shorter time.

When you stop drinking your level of intoxication can change quickly (rise or fall) for up to 90 minutes afterwards because of the time it takes for alcohol to reach the bloodstream after drinking. For example it will take longer for alcohol to reach the bloodstream if you have eaten recently, less time on an empty stomach. There are too many factors that can influence the reading just after drinking, this is the reason why we do not recommend using the breathalyser within 90 minutes of drinking, the ideal time is the morning after.

## Things That Can Affect the Reading of Your Breathalyser

There are some factors that can affect the reading of breathalysers whether it be the AlcoSense breathalyser or a police system. Please note that certain medications and conditions including diabetes can give false positive readings. If you find you have a positive reading when you have not drunk alcohol please contact us for assistance.

### Do:

- Use the breathalyser the morning after a few drinks.
- Read the instructions carefully.
- Phone us if you have any difficulties.
- Only drive when the AlcoSense says 0.00 and you feel 100%. Any amount of alcohol in your system will affect your ability to drive.
- Leave a minimum of 3 minutes between tests. Do not perform more than 10 tests in one hour and no more than 30 tests in 24 hours. More frequent testing can cause a build-up of breath condensation in the breathalyser which can temporarily affect accuracy.

### Don't:

- Blow into the breathalyser within 10 minutes of drinking, high levels of alcohol in the mouth can permanently damage the sensor.
- Do a breath test within 90 minutes of drinking. Blood alcohol levels can rise for up to 90 minutes after drinking (see page 4 for more info). Results taken before this time can rise or fall quickly.
- Use the breathalyser for at least 30 minutes after using mouthwash, breath spray, cough medicine or anything else that may contain alcohol - you may get a false positive reading.
- Use the breathalyser for at least 30 minutes after having a cigarette or eating ripe fruit - you may get a false positive reading.
- Perform a test in dirty air, i.e. a smokey room or pub.
- Use the breathalyser if you are breathing quickly (for example, after running up the stairs) or have been holding your breath. This can falsely lower or raise the reading respectively. Wait until you are breathing normally again.
- Perform two tests sequentially, you may get a variation between readings. This is because it can take a few minutes for the alcohol concentration inside the lungs to stabilise again after the first test.

## Getting Started

This section will take you step by step through the basics of using your new AlcoSense Breathalyser, namely:

1. Unpacking the box
2. Understanding the buttons and screens
3. What the icons mean
4. Installing the batteries
5. First Time Set-Up
6. How to perform a breath test
7. Understanding the results


## Unpacking the Box

In the box you will have the following components:

1. 1 x AlcoSense Breathalyser
2. This Instruction Manual
3. 3 x AAA Batteries
4. 5 x Blow Tubes
5. 1 x micro USB cable for software updates (not charging)



## Understanding the Buttons and Screens

First, take a minute to understand the how to use the screens and buttons on the AlcoSense Excel BEFORE turning it on for the first time. There are three buttons on the front, the ideal place to press each button is shown below by grey circles.

The icons at the bottom of each screen will change according to the function of the button beneath it. For example, in the diagram on the left below, the left hand button accesses the main menu and the middle button will power the unit off.

In the diagram on the right the left hand button will scroll down the menu, the middle button will select the option the white box is highlighting and the right hand button will scroll up the menu.

### What the Icons Mean

Navigation icons are shown at the bottom of the screen, indicating the function of each button.

#### Navigation Icons

Icon	Name	Description
	Power	When you see this icon, press and hold the central button down for 3 seconds to turn the unit off
	Up	Scrolls up a list of options
	Down	Scrolls down a list of options
	OK	Selects an option or progresses you to the next screen
	Main Menu	When you see this icon on the left of the screen, press the left hand button to access the main menu
	Back	This icon will take you back one step
	Repeat / Next	This icon will either allow you to repeat a step or progress to the next step


Message icons are shown on the screen to communicate a message.

#### Message Icons

Icon	Name	Description
	Don't Drive Icon	Displayed when you are close to or over your set limit - do not drive
	Warning Icon	Either signifying you have alcohol in your sample, but are below the limit or alerting you to an error
	Usage Tip	AlcoSense Excel will give you tips to help you perform a good test the next time round.

### Installing the Batteries

Keep the AlcoSense Breathalyser in the closed position and turn the machine round so the back is facing you. Following the diagram below remove the battery compartment cover. Insert the three 1.5v AAA batteries supplied into the slots provided in the battery compartment making sure that the batteries are inserted the right way round as indicated in the battery compartment. Put the battery compartment cover back on again following the diagram below but in reverse.

**IMPORTANT:** When replacing the batteries do not mix old and new batteries. For best accuracy only use a matching set of new good quality branded batteries.

### First Time Setup

1. Press and hold the power button in the middle down until the screen lights up or slide the unit open.
2. Select your language. If you select the wrong language continue with first time setup and see page 19 to see how you can change it back to English
3. **IMPORTANT:** You must set the date correctly to today's date when it asks you. This will be used to calculate when the unit will alert you to recalibrate it. This date cannot be changed once set. Use the up and down buttons to change the numbers, and press the middle button to move to the next item.
4. **Default Drink Drive Limit Settings.** When you first turn your AlcoSense Excel on it is set to measure in %BAC, to alert you at the Scottish drink drive limit of 0.50%BAC (0.22mg/L) and set to the UK Blood Breath Ratio (BBR) of 2300:1. If you wish to change these settings please go to the contents page and find the section called "Changing Limits, Units of Measure and Driving Abroad".
5. **You should now check for any software updates.** Go to the Main Menu > Info and check the "Version" number. Go to [www.alcosense.com/support](http://www.alcosense.com/support) and check if a later version number of software is available for the AlcoSense Excel. If so, follow the instructions on the web page to update your unit. We recommend checking back regularly to get the latest software and keep your unit up to date.

### Turning Off


The unit will enter a power saving mode and dull the screen brightness after 1 minute of inactivity, and turn itself off after a further 2 minutes of inactivity.

If you wish to turn the unit off manually, press and hold the middle key with the power logo on for 3 seconds from any screen.

## How to Perform a Breath Test

AlcoSense is a sensitive measurement device. Performing a breath test less than 10 minutes after drinking risks blowing high levels of alcohol directly from your mouth into the system causing permanent damage to the sensor and causing the system to read inaccurately. Blood alcohol levels can rise for up to 90 minutes after drinking, for best results only use after this period and ideally the morning after. Do not blow directly into the blow hole without a blow tube – this can damage the sensor and significantly reduce the lifetime of your unit.

1. Slide the unit open and insert a blow tube as shown below. Please note that the blow tube sits at an offset angle to make it easier to see the screen while blowing. The tube sits snugly in the recess but is not tight. For the best accuracy and hygiene we recommended using a new blow tube each time.

2. The AlcoSense will display two warnings, please read these carefully. On further tests these warnings can be skipped by pressing the central button twice.

3. You will then see the following screen on the left. You will need to blow for approximately 6 to 8 seconds at a steady pressure (like you're taking a sigh) so make sure you pace your breath for this. When you start blowing into the breathalyser the screen will change to the one in the middle below if the pressure is strong enough, and after 6 to 8 seconds of blowing you will see your result. Keep blowing until you hear a beep.

If you are not blowing hard enough the left hand screen will continue to show. If this is the case, stop blowing and start again with a slightly harder pressure. If you do not blow for long enough or do not blow at a steady pressure the screen on the right will be shown and you must try again.

Insert blow tube

Blow steadily like  
you're taking a sigh

Keep Blowing!

Blow Error  
Blow steadily like  
you're taking a sigh  
OK

### PLEASE NOTE:

Leave a minimum of 3 minutes between tests. Do not perform more than 10 tests in one hour and no more than 30 tests in 24 hours. More frequent testing can cause a build-up of breath condensation in the breathalyser which can temporarily affect accuracy until it evaporates.

If you perform sequential tests you may get varying results. This is due to the alcohol level in your lungs needing a couple of minutes to stabilise after performing the first test, and the fact your body is constantly processing alcohol into and out of your blood stream. As a result UK law allows a 15% difference between two sequential test results.

## Understanding the Results

AlcoSense Excel will alert you not to drive from 20% below your selected drink drive limit. Before using the device, make sure you have set the correct drink drive limit and BBR for your country or region. Please see the relevant sections of this manual which show how to change this and what settings are required.

You can check the current setting in Main Menu > Info, then scroll down to "Limit". Drink drive limits do sometimes change, before making a trip to a different region or country always check the limit with the local authorities and go to the support section on [www.alcosense.com](http://www.alcosense.com) to check for updates.

Description	Screen Example	Alcohol Level
<b>Zero or low level alcohol</b> The reading is below your chosen drink drive limit		Zero to 0.09mg/L (0.20‰BAC).
<b>Caution, Alcohol Present</b> There is enough alcohol in the sample to affect your driving but you are below your chosen limit. We recommend you do not drive		Between 20% and 80% of your set drink drive limit
<b>Close to or Over Limit</b> The sample was very close to, or over the set limit – DO NOT DRIVE		20% below your set drink drive limit or greater

### Sensor Saturation Risk

The unit has received an extremely high sample outside the range expected from a human. This is normally caused by blowing alcohol into the unit directly from the mouth after drinking. Readings taken after this time may be inaccurate until the sensor has recovered. To allow recovery the unit will lock down for 24 hours, and then work as usual afterwards

Above 1.00mg/L (2.30‰BAC)

### IMPORTANT NOTE ON READINGS:

To comply with the pan-European breathalyser standard, EN16280, the AlcoSense Excel adds 20% to the reading shown to the end user. This minimises the risk that the item will under-read even toward the end of the 12 month calibration period. The Excel will also not show a reading below 0.09mg/L (0.20‰BAC) but will instead show the "LOW" screen opposite.

### USING YOUR BREATHALYSER IN A NOISY ENVIRONMENT:


The operation of the breathalyser can be compromised in noisy environments over 95dBA, resulting in a breath test error message being shown. When the unit is slid open to start a test, pressure changes as a result of loud noises can occasionally cause the unit to think breath is entering the device, and then show a blow error. If this occurs please move to a quieter area away from the source of the sound and attempt the test again.

## Changing Settings & Recalling Results

From the home screen, shown below, press the left hand button with the menu icon above it. From here you will find access to previous results, settings and system information.

Once in the menu, the left hand button will scroll down the options, the right hand button will scroll up the options and the centre button will select the option.

### Adjusting the Time & Date

Access the main menu and select "Settings", then select the "Time & Date" option.

### Recalling Last Results

AlcoSense Excel stores your previous 24 test results. When the memory is full the oldest record is deleted to make space for the newest record. To review the previous results navigate to Main Menu > Last Results. Use the up and down key to navigate through the results. The central button will exit back out of this feature.

### Changing the Screen Brightness

From the main screen press the left hand button to access the Main Menu. Navigate to Main Menu > Settings > Brightness.

### Turning the Key Beep On and Off

From the main screen press the left hand button to access the Main Menu. Navigate to Main Menu > Settings > Key Beep.

## Changing the language

You can change the language by going to Main Menu > Settings > Language. If the unit has been set into a different language you can change it back to English by going here:

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## Resetting the Unit

If you ever need to reset the unit back to factory settings navigate to Main Menu > Settings > Advanced > Reset and follow the on screen instructions. The calibration expiry date and your previous test results will not be reset or deleted.

## Service Menu

This option can only be accessed by AlcoSense approved service engineers. Changing settings here will affect the accuracy of your device.

## Changing Limits, Units of Measure and Driving Abroad

### Changing the Drink Drive Limit

This is accessed through Main Menu > Settings > Set Limit. This will change the point at which the AlcoSense Excel will display the don't drive alert. Please note, the don't drive alert is displayed from 20% below your chosen limit. For example, if the limit is set at 0.35mg/L it will display the don't drive alert from 0.28mg/L upwards.

To check the current setting please go to Main Menu > Info and scroll down until you see "Limit". If you wish to change the limit please go to Main Menu > Settings > Set Limit. You can then adjust the limit up or down as you choose.

### Changing the Unit of Measure

You can change the AlcoSense's unit of measurement between blood (BAC) or breath (BrAC) alcohol concentrations. If using a blood alcohol unit of measure we STRONGLY recommend double checking the Blood Breath Ratio (BBR) is set correctly first by referencing the tables at

www.alcosense.co.uk/support. The UK and Irish BBR is 2300:1. If you use the wrong BBR your readings will display incorrectly (see section on Blood Breath Ratio). The choice of units are:

- mg/L BrAC (milligrams of alcohol per litre of breath)
- ‰BAC (per mil Blood Alcohol Concentration)

### Driving Abroad > Changing the Blood Breath Ratio

If you are using your AlcoSense breathalyser abroad you may need to change the Blood Breath Ratio (BBR) setting as well as the drink drive limit. The BBR is the ratio that is used to convert a breath alcohol reading (displayed in mg/L) into a blood alcohol reading (displayed in ‰BAC). Different countries use different BBRs. The UK uses 2300:1 which gives the highest blood value from a breath reading. A more detailed explanation of the Blood Breath Ratio is shown on the page opposite.

A list of the common BBR's is displayed overleaf. If you do not know the BBR in the country you are travelling to, use 2300:1 as this will always give the highest reading, or change the unit of measure to mg/L (breath alcohol) and use the local mg/L limit.

To check the current setting please go to Main Menu > Info and scroll down until you see "Ratio" or "BBR".

To change the BBR go to Main Menu > Settings > Advanced > Blood Breath Ratio.

#### WARNING:

Do not forget to reset the BBR to 2300:1 when you return to the UK. Using a lower BBR in the UK will cause your unit to under-read.

### How can a Breathalyser Display a Blood Alcohol Reading?

Breathalysers measure the alcohol concentration in a persons breath and can convert the breath reading into a blood reading using the Blood Breath Ratio.

#### The Blood Breath Ratio

**IMPORTANT:** The Blood Breath Ratio (BBR) used in drink driving law differs slightly from country to country, but is normally 2000:1, 2100:1, 2200:1 or 2300:1. If the wrong BBR is used in the wrong country, a breathalyser could give a falsely high or low blood alcohol reading. For example, the UK uses a BBR of 2300:1 which means a breath reading of 0.22mg/L = a blood reading of 0.50‰BAC. However, Germany uses a BBR of 2100:1 which means a breath reading of 0.22mg/L = a blood reading of 0.46‰BAC. This means that a motorist who normally drives in the UK, uses the wrong BBR setting abroad could be over-reading elsewhere, or a foreign motorist driving to the UK could under read. This is why it is very important to select the correct BBR in Main Menu > Settings > Advanced > BBR. If you do not know the correct BBR for the country, set it to 2300:1 as this is the option which will give the higher reading.

#### The Relationship Between Different Units of Measure

AlcoSense Excel converts the value for you if you choose a different unit of measure and the correct Blood Breath Ratio is set (see relevant section above). If you wish to convert between units manually the following explains the relationship between different units of measure:

#### Breath Units

0.35 mg/L = 35 µg/100mL

#### Blood Units

0.80‰BAC = 0.080%BAC

To convert a breath reading to a blood reading or vice versa you will need to know the correct Blood Breath Ratio (see relevant section above) for the country you are in. You can convert a breath reading to a blood reading using the following formula:

Reading in mg/L x BBR / 1000 = Reading in ‰BAC

## Common Drink Drive Limits and Blood Breath Ratios

**PLEASE NOTE:** Due to differing Blood Breath Ratios (see page 18–19) used in different countries the same breath reading may equate to a different blood reading in different countries. For example, a breath reading of 0.25mg/L equals a blood alcohol reading of 0.50‰ BAC in Finland, but in Scotland this would equal 0.575‰ BAC. You must ensure you have set the correct Blood Breath Ratio (BBR), if you are not sure of the BBR please set it to 2300:1. See page 18 to see how.w.

Country	Blood Breath Ratio (BBR)	Breath Alcohol	Blood Alcohol
		mg/L Limit*	‰BAC Limit*
Australia	2100:1	0.24	0.50
Australia (Prof. Drivers)	2100:1	0.10	0.20
Austria	2100:1	0.24	0.50
Belgium	2100:1	0.24	0.50
Canada	2100:1	0.24	0.50
Cyprus	2100:1	0.09	0.20
Czech Republic	2100:1	0.00	0.00
Denmark	2000:1	0.25	0.50
Finland	2000:1	0.25	0.50
France (Lower Limit)*	2000:1	0.10	0.20
France (Normal Limit)*	2000:1	0.25	0.50
Germany (Lower Limit)	2100:1	0.00	0.00
Germany (Normal Limit)	2100:1	0.24	0.50
Holland (Lower Limit)	2100:1	0.09	0.20
Holland (Normal Limit)	2100:1	0.24	0.50
Hungary	2100:1	0.00	0.00
Iceland	2100:1	0.14	0.30
Ireland (Rol)	2300:1	0.22	0.50

20

Country	Blood Breath Ratio (BBR)	Breath Alcohol	Blood Alcohol
		mg/L Limit*	‰BAC Limit*
Ireland (Rol) Prof. Drivers	2300:1	0.09	0.20
Italy (Lower Limit)	2100:1	0.00	0.00
Italy (Normal Limit)	2100:1	0.24	0.50
Luxembourg (Lower Limit)	2100:1	0.09	0.20
Luxembourg (Normal Limit)	2100:1	0.24	0.50
Malta	2300:1	0.35	0.80
New Zealand	2100:1	0.24	0.50
New Zealand (Under 20)	2100:1	0.00	0.00
Norway	2000:1	0.09	0.20
Poland	2100:1	0.09	0.20
Russia	2100:1	0.00	0.00
Slovakia	2100:1	0.00	0.00
Spain (Lower Limit)	2000:1	0.15	0.30
Spain (Normal Limit)	2000:1	0.25	0.50
Sweden	2000:1	0.10	0.20
Switzerland (Lower Limit)	2100:1	0.00	0.00
Switzerland (Normal Limit)	2100:1	0.24	0.50
Turkey	2100:1	0.24	0.50
UAE	2300:1	0.00	0.00
UK (England & Wales)	2300:1	0.35	0.80
UK (Northern Ireland)	2300:1	0.35	0.80
UK (Scotland)	2300:1	0.22	0.50
USA (varies by state)	2100:1	0.24 - 0.38	0.50 - 0.80

\*Limits correct as of January 2015, errors and omissions excluded. Always check the limit with the local road authority in the region before driving. Northern Ireland has proposals to reduce their limit, please check with local road authorities before driving. \*\*The AlcoSense Excel cannot be used to comply with French law.

21

## Maintaining your AlcoSense Excel Breathalyser

### Keep it Clean

To keep your AlcoSense in tip top condition we recommend keeping it clean using a soft dry cloth and storing it in a clean, temperate, non dusty environment. Please see the specification section at the back of the manual for further details.

### Be Careful

Take care not to drop the AlcoSense breathalyser as sudden knocks can affect the unit's accuracy. If you have dropped your AlcoSense and feel that it is not reading correctly please call us to arrange to have the unit serviced.

### Keep it up to Date

Regularly check the support section of [www.alcosense.com](http://www.alcosense.com) for software updates. Instructions for updating your unit are available online.

### Servicing

We will happily service your breathalyser to keep it in top condition so that it is reading to the highest levels of accuracy and reliability. Servicing normally takes around 5 working days to complete as after each service your AlcoSense breathalyser will be quality tested before leaving us.

### Please choose your level of service below:

#### 1. Re-Calibration

We will re-calibrate your AlcoSense breathalyser to ensure it continues to read at its most accurate. This service is charged at £24.99 and includes return postage costs to you. Please call 0800 195 0088 (or +44 1628 778885 from outside the UK) to arrange this. We recommend this is done once every 12 months.

#### 2. Out of Warranty Repair

If your breathalyser is out of warranty don't worry, we can repair it for you. You will be charged for the parts and labour required at a reasonable cost.

Please go to [www.alcosense.com](http://www.alcosense.com) to arrange this.

## Technical Support

First, please check the troubleshooting section in this manual and [www.alcosense.com/support](http://www.alcosense.com/support). If you have any queries which are not covered or need assistance using your unit please call our technical support team on 0800 195 0088 (+44 1628 778885 outside the UK) Monday to Friday 9am to 5pm. When calling please make sure you have your breathalyser to hand.

### Re-Calibration

Your AlcoSense breathalyser will require re-calibration after 12 months, after this period there is a risk of obtaining erroneous results. Calibration can only be performed by AlcoSense or an officially authorised service centre. Please contact us to arrange this.

## Troubleshooting

### I have got a positive reading but I haven't drunk any alcohol.

This can be caused by testing within 10 minutes of eating some strong foods, using dental products, chewing gum or smoking. On some occasions it can also be caused by the body releasing acetone/keytones released during the metabolism of acetoacetic acid. This can be increased with the use of some prescribed medications or certain medical conditions, most commonly diabetes. A healthy body can also create these substances occasionally, it is more likely if you are on a low carbohydrate diet such as the Atkins diet or similar. If you experience this problem please call us.

### I'm doing a breath test and I can't get it to work.

This can be because you are blowing too hard, too soft or not for long enough into the mouthpiece. Try blowing like you're taking a sigh. If the screen shows the "insert a blow tube" message still, try increasing the pressure until it says keep blowing. Blow steadily at this pressure until the result is shown. Conversely if you are running out of breath before the result is shown, reduce the pressure of your breath slightly.

### **I've done two breath tests and got different results.**

This is expected for several reasons. The first test is the most accurate, if you perform a second test within a few minutes of the first test the alcohol levels in your lungs may not have stabilised after the first test, causing the second test to give you a different reading. If you have recently drunk alcohol, your blood alcohol level can also change quickly over the course of a few minutes. UK law allows for a 15% difference in sequential readings as a result. This is why you should only use the breathalyser at least 90 minutes after drinking or ideally, the morning after.

Additionally, this is why you cannot use a human to test the accuracy of two different breathalysers back to back, you will also expect slightly different readings for the same reason. It can only be done using approved Breath Alcohol Simulators.

### **My unit says "Saturation Lockdown"**

This is caused by blowing in an abnormally high level of alcohol. It is normally caused by doing a breath test just after sipping an alcoholic drink and there is residual alcohol left in the mouth. This overloads the sensor and can make it read inaccurately for a period of time while the sensor recovers. For your safety the unit stops you from performing further tests until the sensor has recovered. After the period has expired your unit should work as usual again.

### **My unit says "Calibration Expired" or "Calibration Due"**

The 12 month calibration interval is about to or has expired. After this time readings may become unreliable. Please contact us to arrange re-calibration of the unit, our contact details can be found on the back cover of the manual.

### **My calibration expiry date appears to be wrong**

This happens if you do not set the date correctly when you turned the unit on for the first time during initial setup. Please contact us to resolve the issue.

### **I've done a lot of tests in quick succession and I've got a strange result**

This can be caused by a build-up of breath condensation in the breathalyser which sometimes cause temporary inaccuracy until it evaporates, this can affect any breathalyser, including Police breathalysers. To avoid this leave a minimum of 3 minutes between tests and do not perform more than 10 tests in one hour, and no more than 30 tests in 24 hours. Always use a new mouthpiece for each test. If you do experience this problem the condensation should evaporate within 2 hours at normal room temperature and the breathalyser will return to normal operation.

### **I got a LOW Result**

A LOW result means that there is very little or no alcohol in your sample. To comply with the pan-European breathalyser standard, EN16280, the AlcoSense Excel cannot show a reading when a sample of less than 0.09mg/L (0.20‰BAC) is detected. This is remove the likelihood of a false positive result being displayed. Please see the section titled "Understanding the Results" for more information.

### **My breathalyser started to do a test by itself and displayed an error**

This can happen if the breathalyser is used in a noisy environment over 95dBA. When the unit is slid open to start a test, pressure changes as a result of loud noises can occasionally cause the unit to think breath is entering the device, and then show a blow error. If this occurs please move to a quieter area away from the source of the sound and attempt the test again.

### **Any other issues**

Please check the [www.alcosense.com](http://www.alcosense.com) support section for any software or firmware updates, if available please download and update your device. If you still experience problems after this point please contact us.

**For any other queries please call us on 0800 195 0088 from inside the UK, or on +44 1628 778885 from outside the UK.**

## Specification

Reading Range	0.00mg/L to 1.00mg/L, results below 0.04mg/L will display as "LOW" (0.00mg/L)
Accuracy	-0.0%BAC, +0.15%BAC at 0.50%BAC or -0.00mg/L, +0.07mg/L at 0.22mg/L
Operating Temperature	5°C to 40°C
Storage Temperatures	-20°C to 70°C
Storage conditions	Store in a cool dry place out of direct sunlight
Calibration Intervals	12 months
Sensor	64mm <sup>2</sup> Platinum Fuel Cell
Sampling System	Pressure and volume

## Accessories

The following accessories are available for your AlcoSense unit, which are available either from your local dealer or directly from [www.alcosense.com](http://www.alcosense.com)

Accessory	Description	Price*
Protective Carry Case	Keep your AlcoSense breathalyser safe and sound in this beautiful carry case	£12.99
Pack of 20 Spare blow tubes	Pack of 20 blow tubes for breath tests	£4.99
Pack of 25 one-way mouthpieces	Mouthpieces designed for ultimate hygiene so vapour from a previous test cannot enter the mouth of the next user.	£9.99

\*Prices may vary by retailer and are subject to change if bought from AlcoSense directly.

### Note on Environmental Protection

After implementation of the European Directive 2002/96/EU in the national legal system, the following applies: Electrical and electronic devices may not be disposed of with domestic waste. Consumers are obliged by law to return electrical and electronic devices at the end of their service lives to the public collecting points set up for the purpose or point of sale. Details to this are defined by the national law of the respective country. This symbol of the product, packaging or instruction manual indicates that this product is subject to these regulations. By Recycling, re-using the materials or other forms of utilising old devices, you are making an important contribution to protecting our environment.

Although every effort has been made to ensure that the contents of this manual are accurate and reliable Now Group UK Ltd accept no liability for errors or omissions and we reserve the right to change information without notice.

### AlcoSense limited warranty

#### WARRANTOR

This Limited Warranty is granted by and this Limitation of Liability is stipulated for the benefit of AlcoSense, a trade mark of Now Group UK Ltd, at the address displayed on [www.alcosense.co.uk](http://www.alcosense.co.uk)

#### WHAT THIS WARRANTY COVERS

1) Now Group UK Ltd. ("AlcoSense") warrants to you that the Hardware will be free from defects in workmanship and materials under normal use ("Defects") for a period of one (1) year (unless the country of purchase requires a longer period under local legislation) from the date that the Hardware was first purchased by you from an authorised AlcoSense dealer ("Warranty Period"). During the Warranty Period the Hardware will be repaired or replaced at AlcoSense's choice ("Limited Warranty") without charge to you for either parts or labour. This Limited Warranty covers the replacement of the Hardware only. If the Hardware is repaired after the Warranty Period has expired (at your cost), the Warranty Period for the repair will expire six (6) months after the date of repair.

#### WHAT THIS WARRANTY DOES NOT COVER

2) The Limited Warranty does not apply to normal wear and tear, does not apply when the Hardware is opened or repaired by someone not authorised by AlcoSense and does not cover repair or replacement of any Hardware or part thereof damaged by: misuse, moisture, liquids, proximity or exposure to heat and accident, abuse, non-compliance with the instructions supplied with the Hardware, neglect or misapplication. The Limited Warranty does not cover physical damage to the surface of the Hardware or calibration of the Hardware.

#### HOW TO MAKE A WARRANTY CLAIM

3) In order to make a claim of a Defect, you must contact AlcoSense by email or telephone during the Warranty Period via [www.alcosense.com](http://www.alcosense.com) to explain the Defect and to obtain an RMA number (Return Materials Authorisation) if necessary. You must return the Hardware during the Warranty Period, along with an explanation of the Defect, to the address provided to you by AlcoSense. If a defect arises and a valid claim under this Limited Warranty is received by AlcoSense after the first one hundred and eighty (180) days of the Warranty Period, AlcoSense is entitled to charge you for any reasonable shipping and handling costs made in connection with the repair or replacement of the Hardware. You must comply with any other return procedures stipulated by AlcoSense, if any.

#### YOUR LEGAL RIGHTS

4) Some countries may not allow the exclusion or limitation of damages. If any part of this Limited Warranty is held to be invalid or unenforceable, the remainder of the Limited Warranty shall nonetheless remain in full force and effect.  
5) This Limited Warranty is the only express warranty made

to you and is provided in lieu of any other express warranties or similar obligations (if any) created by any advertising, documentation, packaging, or other communications.

6) Except for the Limited Warranty and to the maximum extent permitted by applicable law, AlcoSense and its suppliers provide the Hardware "AS IS AND WITH ALL FAULTS", and hereby disclaim all other warranties and conditions, whether express, implied or statutory, including, but not limited to, any (if any) implied warranties, duties or conditions of satisfactory quality, of fitness for a particular purpose, of reliability or availability, of accuracy or completeness of responses, of results, of workmanlike effort, of lack of viruses, and of reasonable care and skill, all with regard to the Hardware, and the provision of or failure to provide support or other services, information, software, and related content through the Hardware or otherwise arising out of the use of the Hardware. Also, there is no warranty or condition of quiet enjoyment, quiet possession, or non-infringement with regard to the Hardware. This exclusion does not apply to (i) any implied condition as to title and (ii) any implied warranty as to conformity with description.

7) This Limited Warranty does not affect any legal rights under applicable national legislation governing the sale of consumer goods.  
8) This Limited Warranty cannot be transferred to any other person.

#### LIMITATION OF LIABILITY

1) Neither Now Group UK Ltd nor its suppliers shall be liable to you or to any third party for any damages either direct, indirect, incidental, consequential or otherwise (including in each case, but not limited to, damages for the inability to use the equipment, loss of business, loss of profits, business interruption or the like) arising out of the use of or inability to use the Hardware even if Now Group UK Ltd has been advised of the possibility of such damages.

2) Notwithstanding any damages that you might incur for any reason whatsoever (including, without limitation, all damages referenced herein and all direct or general damages in contract or anything else), the entire liability of Now Group UK Ltd and any of its suppliers shall be limited to the amount actually paid by you for the Hardware.

3) Now Group UK Ltd shall not be liable for (i) any fraud on the part of its employees and/or agents; or (ii) any fraudulent misrepresentation on the part of its employees and/or agents.  
4) Notwithstanding the above, neither party's liability for death or personal injury resulting from its own negligence shall be limited.

5) If any of these terms shall become or be declared illegal, invalid or unenforceable for any reason whatsoever such term shall be divisible from these terms and shall be deemed to be deleted from them.

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Patents pending [www.alcosense.com/ip](http://www.alcosense.com/ip) © Now Group UK Ltd
### **Part of Now Group UK Ltd**