PL

EN

BEAUTIFLY

B-Glossy

PRO

series

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yellow, blue and red photon light		
5 working modes	EMS currents	cryotherapy
collagen synthesis function	vibration frequency 8000 / min	
reduces the signs of skin aging		lifts and tightens the skin
narrows pores	new generation PRO facelift	

The B-Glossy PRO needle-free mesotherapy device is a new generation device for anti-aging treatments and many others.

Find out about the advantages of B-Glossy PRO and take a moment to visit your private beauty zone. Take advantage of many possibilities depending on what your skin needs at the moment. Deep intracellular hydration, cryotherapy for swollen tired eyes, and maybe something more advanced? EMS microcurrent will stimulate the skin to raise its oval shape and give it firmness and tension. The device has the Collagen Synthesis function, i.e. a technology that will immediately stimulate the skin to increase collagen production. Collagen relaxes, smoothes wrinkles, tightens the skin, giving a lifting effect in just a dozen minutes. Ultrasonic micro-massage in combination with photon light therapy will deeply penetrate the skin layers, providing it with the maximum amount of active substances. Light therapy will individually affect specific skin problems, such as acne, wrinkles, or a decrease in intracellular hydration, restoring the skin to a naturally healthy, and younger appearance.

Technical data

input:	DC 5V
temperatur:	6-42°C (+/-5°C)
vibration:	8000 (-+ 10) rpm
charging time:	2.5h
charging current:	500mA
power:	max 8.5W
battery capacity:	1000mAh
standby time:	45 days
net weight:	145g
use time:	5min/time or 2 min/time
device size:	175*55*67mm

Introductions

Our complexion consists of at least three million pores. It is known that in order to thoroughly clean it, the liquid and hands alone are not enough. It is best to remove make-up and surface impurities and clean the skin with sonic brushes. From time to time, however, each skin needs a thorough intracellular cleansing. Then the skin will be able to receive and absorb the active ingredients that we apply, and the cosmetics will fulfill their purpose perfectly. Face-adapted sonic micro-vibration technology cleans the skin with 8,000 vibrations per minute. All impurities are carefully brought to the surface and scraped off the face with a firmly attached cotton pad. After such cleansing, the skin is smooth and the pores are cleansed. Ready for further care steps.

The principle of heat cold

The principle of warm and cold Research shows that reaching the skin at 42°C helps to produce HSP protein, repair the skin's nerve cells, and increase skin elasticity to counter the formation of wrinkles. A heat treatment at 42°C helps to relax the pores of the skin for deep cleansing, while a cold treatment with a cold effect to 6°C immediately contracts the pores and blocks nourishment. It alternates between heat and cold care, stimulates the skin, improves blood circulation, and accelerates the metabolism of residual fluids. In addition, it effectively eliminates dark circles under the eyes, smoothes fine lines and puffiness, and also brightens the skin. In combination with micro-vibrations, it has a positive effect on skin firmness and reduces fatigue. The skin becomes radiant and radiant in a natural way.

EMS (Electrical Muscle Stimulation)

EMS contributes to the activation of skin cells by stimulating it to produce collagen, which improves skin elasticity, smoothes wrinkles, and leaves the skin smooth and healthy.

Crio Effect

1. Supports skin firmness, maintains elasticity, prevents allergies and other skin problems.
2. Shrinks pores retain skin moisture and smooth wrinkles.
3. It regenerates muscles, relieves skin fatigue, the feeling of tension, reduces redness, pain, and inhibits the multiplication of acne lesions.

Photon light therapy

1. Red light: 620nm (± 10 nm) wavelength 18 LEDs deeply penetrate the skin 630nm to promote blood circulation, and collagen regeneration improves skin elasticity to keep skin smooth and moisturized. Suitable for flabby and wrinkled skin.
2. Blue light: wavelength 470 (± 10 nm) 18 LEDs deeply penetrate the skin 480nm to kill bacteria, promote blood circulation, can penetrate up to 0.25mm under the bark, destroy the living environment of bacteria, anti-inflammatory, bactericidal, remove acne effectively, acne scars and so on.
3. Yellow photon light: wavelength 600nm (± 10 nm) 18 LEDs deeply penetrate the skin 610nm to improve cell oxygenation. It stimulates lymphatic drainage. Improves the condition of rough skin, used to treat skin redness.

Detailed instruction manual

1. Before first use, please fully charge the device. Charging time is approximately 2.5 hours. To turn the device on, press the "5" button for approx. 3 seconds. By default, the cleaning mode is set to intensity level 1, then the instrument works with vibration and heat. Depending on the individual needs of the skin, the temperature can be adjusted with the "7" button. You can select the photon lights by pressing the "6" button in each mode (5 minutes for mode 1).
2. Press the Mode button to enter the Cool mode (Cryotherapy), the device vibrates and cools. According to individual needs, you can press 7 to adjust the temperature (Crio session lasts 2 minutes. After turning off the fan continues to run for 2 minutes to cool down)
3. Press the power button "5" to enter the EMS electrostimulation mode. The device vibrates and produces a high EMS microcurrent. According to your individual needs, you can adjust the EMS level by pressing the "7" button (default time setting is 5 minutes for mode 3).
4. Press the "5" button to activate the humidifying mode. The device vibrates by generating heat and EMS micro-currents deep into the skin. According to your individual needs, you can adjust the EMS intensity by pressing the "7" key (default time setting is 5 minutes for mode 4).

5

MODE

7

LEVEL

6

LED

5. Press the power button "5" to enter collagen synthesis mode. The device vibrates and cools, emitting EMS microcurrent to the skin. According to your individual needs, you can adjust the EMS intensity by pressing the "7" button (default time setting is 2 minutes for mode 5).

Clean mode

(opens pores and exports dirt from the deep layers of the skin)

Before using B-Glossy PRO, wash your makeup thoroughly and clean your skin of all impurities. You can use a sonic toothbrush for this. To use the Clean mode, a cotton swab will be required. Remove the transparent rim and place a cotton swab soaked in micellar or another makeup remover in the center of the head. Then put a transparent hoop on it to fix the swab. Start the device by entering the purge mode. Your skin will be heated to open the pores and transport dirt to the outside, leaving it on the swab. Your pores will remain cleansed and your skin will be ready for the next steps of care.

Cool mode

(reduces and closes pores, tightens the skin)

Cool mode or Cryotherapy effectively closes the pores and maintains hydration in the skin. The temperature from 6-42°C alternately makes the skin firmer in a natural way, it becomes smoother and more compact.

EMS tightening mode

(stimulates the muscle cells and tightens the skin without wrinkles)

Before starting the EMS mode, apply a serum or an ampoule suitable for your skin's needs to your face. EMS micro-currents stimulate the muscles, make the skin elastic, and energized. They remove crow's feet and wrinkles on the surface of the skin and restore its youthful appearance, firmness, and smoothness.

Moisturizing mode

(Facilitates the absorption of nutrients into the skin, lifts and tightens)

After cleansing, apply a skincare product such as a serum, ampoule, or essence to your face. Activate the moisturizing mode to transport active ingredients deep into the skin. In each mode, you can activate red, blue, or yellow photon light depending on the skin's needs.

Collagen synthesis mode

(Stop nourishment in the skin, activate collagen synthesis)

After absorbing skincare products, turn on collagen synthesis mode, cool down to shrink pores, and retain skin moisturizing and care. EMS microcurrents will stimulate the skin cells to synthesize collagen and tighten the muscles, giving the effect of skin tension and lifting. In each mode, you can choose red, blue, or yellow photon light according to your skin needs.

Safety instructions

1. Read all safety instruction before using.
2. Follow the operating instructions to use this device.
3. Do not pull down or remake the device, please consult the distributor for repairing.
4. Only use the device for its intended purpose as described in this manual.
5. Do not store or use the device in the high temperature and humidity place such as under sunshine or near fire or in the bathroom, etc. To avoid fire or electric shock. Do not use or place it in or near the high temperature appliance.
6. Do not use and place the device in the humid condition, in case of any short circuit, electrical shock or fire.
7. Keep away from children and physical disabled people allow to use, close supervision is necessary when this device is operated by user under 18 years of age.
8. When operate cool function, the power consumption is relatively large, the fan will have a delay in heat dissipation. It needs to be placed for about 10 minutes for next time use after use cool function.
9. Do not leave the device connected to a power supply for a long time.
10. Do not touch it with wet hands; it may result in electrocution or short-circuit.
11. Do not leave in a high humid places; it may result in electrocution, short-circuit or damage of insulation.
12. Keep out of reach of children.
13. The device may not be used by children.
14. It is forbidden to use the device in places where flammable materials are stored.

15. Do not clean with caustic substances.
16. Use only the adapter provided with the device.
17. Before connecting make sure that the device is dry.
18. Do not use the device on sunburnt skin.
19. Turn off the device in case of any undesired side effects.

Warning

1. Please stop using when the product becomes excessively hot to avoid burns.
2. Do not drop or hit the device.
3. Do not use the device if the power cord is damaged.

Contraindications

- pregnancy and lactation
- metal implants, pacemaker
- metal items (e.g. jewellery)
- cancer
- fever
- medical condition, weakness
- skin breakdown
- use of fillers and botox
- viral, bacterial and fungal skin diseases
- herpes
- epilepsy
- migraines
- thyroid disease
- acute rosacea and phlegmonous acne

Indications

- loss of skin tightness
- grey and sallow skin
- wrinkles
- discolourations
- dry, dehydrated skin

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The crossed-out wheeled bin symbol indicates that in the European Union after the end of the product's useful life it must be disposed of in a separate, dedicated collection point. Do not dispose of these products together with unsorted general household waste

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Plac Bankowy 2, 00-095 Warszawa

English

Hereby, Beauty Factor Sp. z o. o. declares that the device equipment Beautify B-Glossy PRO is in compliance with Directive 2014/30/UE.

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TERMS OF WARRANTY

1. The product warranty is valid for a period of 24 months.
2. The warranty covers the period from the date of receipt of the device from Beauty Factor Sp. z o. o. and covers defects resulting from causes inherent in the received device.
3. The device will be repaired as soon as possible, not exceeding 21 working days from the date of acceptance of the device for repair by an Authorized Service Center.
4. The customer has the right to apply for replacement of equipment free from defects, if the Authorized Service finds that the removal of the defect is impossible or requires excessive costs.
5. If only part of the product is defective and can be detached from the part of the product working properly, the Customer's rights resulting from this warranty limit only to the defective part of the product.

In order to benefit from the warranty conditions, please send a complaint notification together with the warranty card to the e-mail address serwis@rgmedia.pl. The exact warranty conditions are under the address:

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NR DOWODU ZAKUPU
NUMBER OF EVIDENCE OF PURCHASE

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