optimum
med♡

Poduszka do masażu
HZ-7020

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Waga: 1,16 kg 

Dear customers,

Thank you very much for choosing one of our products. Please read these operating instructions carefully and keep them for future reference, share them with other users and follow the information contained therein. Always follow all instructions and basic precautions when using the product.

The device is for private use; it is not designed for medical or industrial use.

WARNINGS AND SAFETY NOTES

- This device may be used by children over the age of 8 and by people with reduced physical, sensory or mental skills or a lack of experience or knowledge, provided that they are supervised or have been instructed on how to use the device safely, and are fully aware of the consequent risks of use.
- Children must not play with the device.
- Cleaning and user maintenance must not be performed by children unless supervised.
- **WARNING!** Keep the packaging away from children. Danger of suffocation.
- During operation, the outer surface becomes hot. People who are sensitive to heat must use the device carefully.
- Never use this appliance if it has a damaged cord or plug, is not working properly, has been dropped or damaged, or has been dropped into water. Return the device to the service center for inspection and repair.
- Under no circumstances should you open or open yourself repair the device, as this may prevent it from working properly action. Failure to comply with the above rules will result in loss warranty.
- The device should only be used with the supplied power supply and using the current characteristics marked on the power supply.
- Do not use when the device is visibly damaged,
- Do not use during a thunderstorm.
- After completing the massage and before each cleaning, turn off the device and remove the plug from the socket.
- Do not pull on the power cord or device to disconnect the power adapter from the outlet.
- Never hold or carry the device by the power cord.
- Do not squeeze, kink or twist the cable.
- Make sure that the massage device, switch, power supply and cable do not come into contact with water, steam or other liquids. Only use the device in dry interiors (e.g. never in a bathtub or sauna).
- Never reach for the device that has fallen into water.
- Do not use under a cover such as a blanket, pillow, ...
- Do not use near gasoline or other flammable materials.
- Repairs to electrical devices may only be carried out by qualified personnel. Incorrectly carried out repairs can seriously endanger the user.

- Protect the device from high temperatures.
- Use this device only for the purpose described in this manual.

APPLICATION

The massage device has rotating massage heads that mimic the squeezing and kneading movements of a real Shiatsu massage. These heads move synchronously and opposite - to the left and right. The massage pillow has the ability to turn on or off the heating function. The heat emitted by the device is designed to warm and additionally relax the mass of muscles.

This device is designed to massage individual parts of the human body:

- relieves pain (e.g. neck, back, spine)
- relaxes tense and contracted muscles
- accelerates muscle regeneration
- relaxes and reduces fatigue
- makes it easier to fall asleep
- soothes and calms the nerves
- helps in the fight against cellulite

If you are not sure if the massage device is suitable for you, please consult your doctor.

WARNING

Do not use the device

- in the event of lesions or injuries in the place to be massaged (e.g. disc prolapse, open wounds),

- massage in the area of the heart, when you have a pacemaker, and to other parts of the body after consulting a doctor,
- for swollen, burned parts of the body, injured parts of the body or skin, or that show inflammation,
- in lacerations, capillary wounds, varicose veins, acne and other skin diseases,
- during pregnancy,
- on the face (eyes), larynx or other particularly sensitive parts of the body,
- during sleep.

Before using the device, consult a doctor:

- if you have a severe disease, cancer or have undergone surgery on the upper body,
- with pacemakers, implants or other assistive devices,
- with thrombosis,
- with pains of unknown cause.

APPLIANCE DESCRIPTION

The set includes AC and 12 V car adapter

DEVICE INSTRUCTION MANUAL

The pillow can be placed under the neck or placed on the backrest of the chair. We can place it under the sacral spine in a lying position or lean against it with the loins or thighs while sitting. We can also put our hands or calves on the pillow.

1. Connect the product to a power source.
2. Press the ON / OFF button to turn on the device – the heating function is active.
3. Press the power button again to change the rotation direction of the massage heads.
4. Press the power button again to turn off the heating function.
5. To turn off the device, press the power button again.

After use, disconnect the product from the power supply.

The massage device should be used for a maximum of 10-15 minutes / 1 cycle.

The device will turn off automatically after 15 minutes.

CLEANING AND MAINTENANCE

- Before cleaning, make sure that the product is turned off and disconnected from the power supply.
- Please clean with a dry cloth.
- Keep the product clean and away from humid environment.
- Do not use cleaning agents that contain solvents.

MAINS OPERATION:

Voltage 100-240 V~; 50/60 Hz

DC 12 V --- 2 A

Rating power: 24 W

ECOLOGY – ENVIRONMENTAL PROTECTION

The “crossed-out trash bin” symbol on electrical equipment or packaging indicates that the device cannot be treated as general household waste and should not be disposed of into containers for such waste.

Obsolete or broken-down electrical device should be delivered to special designated collection points, organized by local public administration, whose purpose is to collect recyclable electrical equipment.

This way, each household helps reduce possible negative effects influencing natural environment and allows the reclaiming of materials which the product is made of.

Weight: 1.16kg