## ENGLISH

### TECHNICAL SPECIFICATIONS

- Input: DC 5V / 2 A
- Battery: 3.7 V / 1800 mAh / 15 W
- Functions: suction and heat
- Magnet therapy
- Accessories: tweezers, spare filters and caps

### INSTRUCTIONS

- Charge the Suction & Heat Anti-Cellulite Massager InnovaGoods® in a USB source.
- Clean the skin and apply plenty of body oil to it.
- Before each use, insert a filter with the help of the tweezers and place the cap (fig.1). Place the second filter (fig.2) to prevent the lotion from penetrating the massager.
- Press the power button; select the desired heat level with the left arrow and the desired suction level with the right arrow.
- **IMPORTANT NOTE:** start the massage with the lowest suction level and gradually increase it, little by little, according to your needs; the magnet therapy is continuous.
- Glide the Suction & Heat Anti-Cellulite Massager InnovaGoods® over your skin using slow and continuous movement from top to bottom on legs, glutes, waist, abdomen and arms.
- Once the massage is completed, switch off the device; it automatically switches off after 25 minutes of use.
- Remove the filters (1 and 2) and clean them with warm water after each use; leave them to air dry before using them again; when they deteriorate, throw them away and use new replacement ones.

### WARNINGS

Do not use on varicose veins. Do not use during pregnancy or menstruation. Do not use it for more than 25 minutes per session. Leave the skin to rest for a few days before the next session.

Fully charge the battery before its first use. Make sure that the voltage of the power source corresponds to that of the appliance. Do not pull on the cable and ensure it does not tangle. Keep the cable away from hot and sharp surfaces. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance or cable in liquid. Do not handle it with wet hands. Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the appliance unsupervised when it is turned on. If the appliance, the cable or the pieces looked damaged or do not work correctly, do not use them or attempt to repair them; they should only be repaired by a qualified technician. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling. Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help

prevent harmful consequences to personal and environmental safety.

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ALL RIGHTS RESERVED. MADE IN P.R.C.