BEAUTIFLY

ComfortZen MoxiTherapy System

FitRecovery

series

Moxibustion Thermal Therapy Seat

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- www.beautify.eu

- www.youtube.com/beautifyeu

| | |
|-------------------|---------------------------|
| Voltage: | 7,4 V |
| Rated Power: | 35 W |
| Battery Capacity: | 6000 mAh |
| Charging Time: | approximately 180 minutes |
| Battery Life: | approximately 90 minutes |

Safety

Before using the product, please read this manual and follow the instructions while using the device.

Display Description [Page 2 - f1]

1. Timer
2. Vibration
3. Heating
4. Moxibustion
5. Red Light / On/Off Button
6. Pillow
7. Moxibustion Area
8. Display

How to Use

1. Place the product on a chair or flat surface with the cover open.
2. Sit on the device.
3. Press and hold the On/Off button for 2 seconds to turn on the device. To turn it off again, press and hold the button for 2 seconds.
4. Choose the mode you want to use:
 - Timer Function - In the device's on state, press the button shortly to activate the function. Short pressing it again adds 15 minutes to the existing time, with a maximum limit of 60 minutes. Long pressing the button resets the time to 15 minutes.
 - Vibration Function - In the device's on state, press the vibration button shortly to activate the function. Short pressing it again allows you to adjust one of the 5 levels of vibration. Long pressing the "vibration" button turns off the function.
 - Heating Function - In the device's on state, short pressing the heating button activates the heating function. Short pressing it again allows you to adjust one of the 5 available heat levels. Long pressing the button turns off the heating function.
 - Moxibustion Function - Lift the cover [Page 2 - f1, point 7] and place moxa. Close the protective cover. In the device's on state, press the button shortly to activate the function. Short pressing it again allows you to adjust one of the 5 levels. Long pressing the button turns off the function.
 - Red Light Function - In the device's on state, press the button shortly to activate the red light. Short pressing it again turns it off.
5. After use, fold the foldable cover back to its original state.

Note: Maximum weight capacity is 120 kg.

Cleaning and Maintenance:

Before cleaning, turn off the power and wait for the device to cool completely. The product is not suitable for washing and ironing. To clean the soiled surface, use a soft cloth soaked in a neutral detergent. Do not use abrasive liquids, abrasive agents, brushes, furniture polish, or paint thinners. Do not immerse in any cleaning liquid to avoid short-circuiting and electrical shock.

Hazards and Warnings:

The product should be used with clothing to avoid potential allergies or irritation. Do not insert pins or other sharp objects that may damage the device. Store away from water, heat, extreme cold, in a dry, well-ventilated place. Do not place heavy objects on it to prevent damage. In case of malfunction, contact technical personnel. Do not use replacement parts that are not approved by the manufacturer. Non-professional service personnel should not disassemble and replace parts of the device without permission. The product is not a toy - be cautious of potential burns and keep it out of the reach of children. The device cannot be turned on while charging. If not used for an extended period, it is recommended to charge it every month. The device is not suitable for pregnant women and individuals with a pacemaker.

For more information about the company and products [Page 2 - f2]

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after 13.08.2005.

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[EN] The manufacturer declares under its sole responsibility that the ComfortZen MoxiTherapy System product is in compliance with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, Directive 2014/53 of the European Parliament and of the Council /EU of 16 April 2014 on the harmonization of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances. The full content of the declaration of conformity is available at: www.beautify.eu/deklaracje

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